# User Guide 1

## Safety manual for hearing aids

for US market

# Content

|                                             |    |
|---------------------------------------------|----|
| Safety information                          | 3  |
| Intended purpose / (use)                    | 9  |
| Explanation of symbols                      | 10 |
| General warnings                            | 12 |
| For certain hearing aid types               | 22 |
| For children under the age of 3 years       | 26 |
| For certain battery types                   | 27 |
| When using remote control apps              | 28 |
| For Hearing Care Professionals              | 29 |
| Important information                       | 31 |
| Operating, transport and storage conditions | 31 |
| Disposal information                        | 32 |
| Conformance information                     | 32 |
| For US market only                          | 40 |
| For prospective hearing aid users           | 40 |
| For hearing aid users                       | 42 |
| Tinnitus feature                            | 46 |
| Cybersecurity                               | 49 |
| For children with hearing loss              | 50 |
| Warranty and service                        | 51 |
| Support                                     | 52 |

# Safety information

This safety manual provides safety information and other important information about your hearing aids. It covers several hearing aid types and optional features.

- ▶ Refer to User Guide 2 (Instructions for Use) for your hearing aids to check the hearing aid type and the activated features.

## **WARNING**

**People younger than 18 should go to a doctor before using this.**

People younger than 18 years old need specialized care, and using this without a medical evaluation may worsen impairment or disability. A hearing aid user who is younger than 18 should have a recent medical evaluation from a doctor, preferably an ear-nose-throat doctor (an ENT). Before using this, a doctor should determine that the use of a hearing aid is appropriate.

## **WARNING**

### **to Hearing Aid Dispensers:**

You should advise a prospective hearing aid user to consult promptly with a doctor, preferably an ear specialist such as an ENT, before dispensing a hearing aid if you determine through inquiry, actual observation, or review of any other available information concerning the prospective user, that the prospective user has any of the following conditions:

- ▶ Visible deformity of the ear, either congenital or traumatic
- ▶ Fluid, pus, or blood coming out of the ear within the previous 6 months
- ▶ Pain or discomfort in the ear
- ▶ History of excessive ear wax or suspicion that something is in the ear canal
- ▶ Dizziness, either recent or long-standing
- ▶ Sudden, quickly worsening, or fluctuating hearing loss within the previous 6 months
- ▶ Hearing loss or ringing (tinnitus) only in one ear or a noticeable difference in hearing between ears
- ▶ Audiometric air-bone gap equal to or greater than 15 dB at 500 Hz, 1000 Hz, and 2000 Hz

**WARNING****to Hearing Aid Dispenser, Outputs over 132 dB SPL:**

You should exercise special care in selecting and fitting a hearing aid with a maximum output that exceeds 132 dB SPL because it may impair the remaining hearing of the hearing aid user.

**CAUTION****This is not hearing protection.**

You should remove this device if you experience overly loud sounds, whether short or long-lasting. If you are in a loud place, you should use the right kind of hearing protection instead of wearing this device. In general, if you would use ear plugs in a loud place, you should remove this device and use ear plugs.

**CAUTION****The sound output should not be uncomfortable or painful.**

You should turn down the volume or remove the device if the sound output is uncomfortably loud or painful. If you consistently need to turn the volume down, you may need to further adjust your device.

## **CAUTION**

**You might need medical help if a piece gets stuck in your ear.**

If any part of your hearing aid, like the eartip, gets stuck in your ear, and you can't easily remove it with your fingers, get medical help as soon as you can. You should not try to use tweezers or cotton swabs because they can push the part farther into your ear, injuring your eardrum or ear canal, possibly seriously.

## **NOTICE**

**What you might expect when you start using a hearing aid.**

A hearing aid can benefit many people with hearing loss. However, you should know it will not restore normal hearing, and you may still have some difficulty hearing over noise. Further, a hearing aid will not prevent or improve a medical condition that causes hearing loss.

People who start using hearing aids sometimes need a few weeks to get used to them. Similarly, many people find that training or counseling can help them get more out of their devices.

If you have hearing loss in both ears, you might get more out of using hearing aids in both, especially in situations that make you tired from listening - for example, noisy environments.

## **NOTICE**

### **Tell FDA about injuries, malfunction, or other adverse events.**

To report a problem involving your hearing aid, you should submit information to FDA as soon as possible after the problem. FDA calls them “adverse events,” and they might include: skin irritation in your ear, injury from the device (like cuts or scratches, or burns from an overheated batter), pieces of the device getting stuck in your ear, suddenly worsening hearing loss from using the device, etc.

Instructions for reporting are available at <https://www.fda.gov/Safety/MedWatch>, or call 1-800-FDA-1088. You can also download a form to mail to FDA.

## **NOTICE**

### **Hearing loss in people younger than 18.**

- ▶ People younger than 18 should see a doctor first, preferably an ear-nose-throat doctor (an ENT), because they may have different needs than adults.
- ▶ The doctor will identify and treat medical conditions as appropriate.
- ▶ The doctor may refer the person to an audiologist for a separate test, a hearing aid evaluation.
- ▶ The hearing aid evaluation will help the audiologist select and fit the appropriate hearing aid.

A person who is younger than 18 years old with hearing loss should have a medical evaluation by a doctor, preferably an ENT, before buying a hearing aid.

The purpose of a medical evaluation is to identify and treat medical conditions that may affect hearing but that a hearing aid won't treat on its own.

Following the medical evaluation and if appropriate, the doctor will provide a written statement that the hearing loss has been medically evaluated and the person is a candidate for a hearing aid.

---

The doctor may refer the person to an audiologist for a hearing aid evaluation, which is different from the medical evaluation and is intended to identify the appropriate hearing aid.

The audiologist will conduct a hearing aid evaluation to assess the person's ability to hear with and without a hearing aid. This will enable the audiologist to select and fit a hearing aid for the person's individual needs. An audiologist can also provide evaluation and rehabilitation since, for people younger than 18, hearing loss may cause problems in language development and educational and social growth. An audiologist is qualified by training and experience to assist in the evaluation and rehabilitation of hearing loss in people younger than 18.

## **Intended purpose / (use)**

Hearing aids are intended to improve the hearing of hearing impaired persons. Diagnosis and prescription of a hearing aid must be performed by hearing health specialists, e.g. acousticians, audiologists or ENT doctors.

Use the hearing aids and accessories only as described in the respective user guides.

## Explanation of symbols

### Symbols used in this document

Points out a situation that could lead to serious, moderate, or minor injuries.

Indicates possible property damage.

Information and warnings for children with hearing aids.

### Symbols on the device or packaging

CE compliance label, confirms compliance with certain European Directives, refer to section "Conformance information".

EMC and radio communications compliance label Australia, refer to section "Conformance information".

Optional symbol for custom models with wireless functionality.

Indicates the legal manufacturer of the device.

Do not dispose of the device with general domestic waste. Read more in section "Disposal information".

Packaging can be recycled.

Read and follow the instructions in the user guide.

## Symbols on the device or packaging

According to IEC 60601-1, the applied parts of this device are classified as type B.

Keep dry during transportation.

Temperature during transportation. Read more in section "Operating, transport and storage conditions".

Relative humidity during transportation. Read more in section "Operating, transport and storage conditions".

Atmospheric pressure during transportation. Read more in section "Operating, transport and storage conditions".

Medical device

Compatible with Qi-certified wireless chargers.

Bluetooth® low energy technology.

## General warnings

### **WARNING**

Risk of impairing the residual hearing of the user.

- ▶ Use only hearing aids that have been fitted especially for your needs.

### **WARNING**

Risk of injury!

- ▶ Do not use obviously damaged devices and return them to point of sale.

### **WARNING**

Note that any unauthorized changes to the product may cause damage to the product or cause injury.

- ▶ Use only approved parts and accessories. Ask your Hearing Care Professional for support.
- ▶ Do not modify this product without authorization from the manufacturer.
- ▶ If this product is modified, appropriate inspection and testing must be conducted to ensure safe use of the product.

### **WARNING**

Your hearing aids may reduce certain background sounds, potentially also traffic or warning signals.

### **WARNING**

Risk of explosion!

- ▶ Do not use your hearing aids in explosive atmospheres (e. g. mining areas).

**WARNING**

This product can expose you to nickel, which is known to the State of California to cause cancer. For more information, go to [www.P65Warnings.ca.gov](http://www.P65Warnings.ca.gov).

**WARNING**

Choking hazard!

Your hearing aids contain small parts which can be swallowed.

- ▶ Keep hearing aids, batteries and accessories out of reach of children and mentally disabled persons.
- ▶ If parts have been swallowed consult a physician or hospital immediately.

## **WARNING**

Keep batteries out of reach of children.

- ▶ For some types of button batteries, severe or fatal injuries can occur in 2 hours or less, if it is swallowed or placed inside any part of the body.
- ▶ If you suspect your child has swallowed or inserted a button battery, call the 24-hour National Battery Ingestion Hotline immediately at 1-800-498-8666 for fast, expert advice.
- ▶ Examine devices and make sure the battery compartment is correctly secured, e.g. that the screw or other mechanical fastener is tightened. Do not use if compartment is not secure.
- ▶ Dispose of used button batteries immediately and safely. Used/discharged batteries can still be dangerous.
- ▶ Tell others about the risk associated with button batteries and how to keep their children safe.

If you have hearing aids that are intended for the fitting of children under the age of 3 years or persons with a developmental age of under 3 years, refer also to section "For children under the age of 3 years".

**NOTICE**

- ▶ Protect your hearing aids from high humidity\*. Do not wear them in the shower\*\* or when you apply make-up, perfume, aftershave, hairspray or suntan lotion.

**NOTICE**

- ▶ Protect your hearing aids from extreme heat. Do not store them in direct sunlight.

**NOTICE**

- ▶ Do not dry your hearing aids in the microwave oven.

**NOTICE**

Different types of strong radiation, e.g. during X-ray or MRI head examinations, may damage hearing aids.

- ▶ Do not wear the hearing aids during these or similar procedures like electrosurgery (diathermy, electrocautery, HF-surgery).

Weaker radiation, e.g. from radio equipment or airport security, does not damage the hearing aids.

\* Statement "high humidity" is obsolete for water resistant devices!

\*\* Statement "Do not wear them in the shower" is obsolete for water resistant devices!

---

**NOTICE**

Your devices comply with international standards. However, it cannot be guaranteed that all products on the market work interference-free, for example some induction cookers, electromagnetic security systems (e.g. anti-theft systems and metal detectors) cellular telephones, RFID, room surveillance systems may cause audible interference. If you encounter unwanted noise in the vicinity of such device, you should increase the distance between the hearing aid and that device.

---

**NOTICE**

In case of a serious incident, report the occurrence to the manufacturer of the device.

---

**NOTICE**

When not in use, please store your hearing aids in the supplied jewel case or charging case (if applicable) and keep out of the reach of children.

## Contraindications

---

### **WARNING**

In case of a known contact allergy, consult the hearing aid manufacturer regarding the presence of the allergen before using the hearing aid. If an allergic reaction develops after wearing the hearing aid, immediately discontinue use and consult a licensed physician.

---

### **WARNING**

Consult a Hearing Care Professional if you experience any unusual side effects like skin irritation, excessive accumulation of ear wax, dizziness, change in your hearing, or if you think there may be a foreign object in your ear canal.

## For hearing aids with wireless functionality

In some countries restrictions for the usage of wireless equipment exist.

- ▶ Refer to local authorities for further information.

### **WARNING**

Risk of affecting electronic equipment!

- ▶ In areas where the use of electronics or wireless devices are restricted, verify if your device has to be turned off.

### **NOTICE**

Your hearing aids are designed to comply with international standards on electromagnetic compatibility but interference with nearby electronic devices could occur. In this case, move away from the source of interference.

### **WARNING**

Interference with implantable devices

- ▶ The device can be used with all electronic implants that comply with the ANSI/AAMI/ISO 14117:2012 Electromagnetic compatibility standard for human implants.
- ▶ In case your implant does not comply with the ANSI/AAMI/ISO 14117:2012 standard, consult the manufacturer of your implantable device for information about the risk of disturbance.

## For hearing aids with a tinnitus feature

Your hearing aid may be equipped with a tinnitus feature (e.g. with a tinnitus therapy signal, a tinnitus noiser, or another special function). The use of the tinnitus feature should be only on the advice and in consultation with your Hearing Care Professional.

### **WARNING**

Risk of further impairment to the user's hearing health.

The volume of the tinnitus feature can be set to a level which could lead to permanent hearing damage when used for a prolonged period of time.

- The tinnitus feature should never be used at uncomfortable levels.

## For hearing aids with AutoPhone magnet

### **WARNING**

Risk of affecting life support systems!

- ▶ Use a magnet only when it is a safe distance away from life support systems, such as pacemakers or magnetic valves. For example, the safe distance between pacemaker and magnet should be at least 10 cm (4 inches).

### **NOTICE**

Magnets can disturb electrical devices and delete stored data.

- ▶ Keep magnets away from computers, monitors, television sets, storage media and other electronic equipment/devices.

### **NOTICE**

In close proximity, AutoPhone magnets can damage receiver units of RIC hearing aids.

- ▶ Keep a minimum distance of 2 cm (0.8 inches) between the RIC receiver and any magnet, including the AutoPhone magnet.

For example, do not store your RIC hearing aids and a smartphone with attached AutoPhone magnet together in a small pocket or box.

## For hearing aids with a magnet in the battery compartment

### **WARNING**

Risk of interference with active and non-active implants!

There is a magnet in the battery compartment, therefore if you wear an active or a non-active implant, e.g. a brain implant:

- ▶ Prior to use, have the electromagnetic compatibility verified.  
Consult the physician that implanted the device before using your hearing aids.
- ▶ Keep a safe distance of about 4 cm (1.6 inches) between the implant and hearing aids.

## For certain hearing aid types

### BTE models

Hearing aids that are worn behind the ear and that have a tube.

#### **CAUTION**

Risk of injury!

- ▶ Always wear the tube with an ear piece.
- ▶ Make sure that the ear piece is completely attached.

#### **CAUTION**

Risk of injury!

- ▶ Insert the ear piece carefully and not too deeply into the ear.

#### **CAUTION**

Risk of injury!

- ▶ In very rare cases the ear piece could remain in your ear when removing the hearing aid. If this happens, have the ear piece removed by a medical professional.

## RIC models

The receiver is placed within the ear canal and connected to the hearing aid via a receiver cable.

### CAUTION

Risk of injury!

- ▶ Always wear the receiver cable with an ear piece.
- ▶ Make sure that the ear piece is completely attached.

### CAUTION

Risk of injury!

- ▶ Insert the ear piece carefully and not too deeply into the ear.

### CAUTION

Risk of injury!

- ▶ In very rare cases the ear piece could remain in your ear when removing the hearing aid. If this happens, have the ear piece removed by a medical professional.

### NOTICE

- ▶ Do not pull the receiver connection as this could damage your hearing aids.

## Instant fit models

The hearing aid is placed within the ear canal and kept in place with an ear piece.

### CAUTION

Risk of injury!

- ▶ Always wear the hearing aid with an ear piece.
- ▶ Make sure that the ear piece is completely attached.

### CAUTION

Risk of injury!

- ▶ Insert the hearing aid carefully and not too deeply into the ear.

### CAUTION

Risk of injury!

- ▶ In very rare cases the ear piece could remain in your ear when removing the hearing aid. If this happens, have the ear piece removed by a medical professional.

## Custom, IIC, CIC and ITC models

The hearing aid is placed within the ear canal.

### **CAUTION**

Risk of injury!

- ▶ Insert the hearing aid carefully and not too deeply into the ear.

## Other models

For all other models, no model-specific safety information applies.

## For children under the age of 3 years

There are special hearing aids for the fitting of children under the age of 3 years or persons with a developmental age of under 3 years. Ask your Hearing Care Professional for further information.

### **WARNING**

Choking hazard!

Your hearing aids contain small parts which can be swallowed.

- ▶ Ensure adequate supervision if infants, small children or mentally disabled persons need to wear hearing aids.
- ▶ Check the completeness of the hearing aids regularly.
- ▶ Ensure that your child or the mentally disabled person does not detach the hearing aid from the earmold.
- ▶ Consult your Hearing Care Professional if the housing is deformed.
- ▶ Keep the battery compartment locked. Verify the proper function of the locking mechanism.
- ▶ Keep batteries and accessories out of children's or mentally disabled person's reach.
- ▶ If swallowed consult a physician or a hospital immediately.

## For certain battery types

### For hearing aids with built-in power cells (lithium-ion rechargeable batteries)

Read the safety information on power cells in the hearing aid's user guide.

### For hearing aids with replaceable batteries

#### NOTICE

Only use zinc-air batteries or nickel-metal hydride (NiMH) rechargeable batteries.

Do not use e.g. silver-zinc or lithium-ion rechargeable batteries.

#### NOTICE

Leaking batteries damage the hearing aids.

- ▶ Turn the hearing aids off when not in use to preserve the battery.
- ▶ Remove batteries when the hearing aids are not in use for a prolonged period of time.

The following table helps you to find the correct zinc-air battery. The suitable battery size for your hearing aids is listed in the user guide.

| Battery size | IEC code | ANSI code | Typical color code |
|--------------|----------|-----------|--------------------|
| 10           | PR70     | 7005ZD    | yellow             |
| 13           | PR48     | 7000ZD    | orange             |
| 312          | PR41     | 7002ZD    | brown              |
| 675          | PR44     | 7003ZD    | blue               |

## When using remote control apps

When using an app for controlling hearing aids:

### **WARNING**

Risk of hearing damage!

The device with the app for controlling hearing aids generates short control signals which may be audible. If the device running the app has a very high audio output there is the risk of hearing damage.

While using the app:

- ▶ Do not hold the loudspeaker of the device to your ears or the ears of others.
- ▶ Do not use the device with headphones, headsets or other audio playback devices.

## For Hearing Care Professionals

---

### **WARNING**

For hearing aids with an output sound pressure level of 132 dB SPL or more:

Risk of impairing the residual hearing of the user.

- ▶ Take special care when fitting this hearing aid.

### **Tinnitus feature**

The target population is primarily the adult population over 21 years of age. The patient may have some control of the level or volume of the signal and the patient should discuss this adjustment as well as his or her comfort level and sound of the signal with their Hearing Care Professional.

### **WARNING**

Risk of further impairment to the user's hearing health.

The volume of the tinnitus feature can be set to a level which could lead to permanent hearing damage when used for a prolonged period of time.

- ▶ Should the tinnitus feature be set to such a level in the hearing aid, advise the user of the maximum amount of time per day he or she should use the tinnitus feature.

For example, occupational safety guidelines restrict continuous noise exposure of 80 dBA SPL to 8 hours per day.

- ▶ The tinnitus feature should never be used at uncomfortable levels.

For hearing aids that are not specially designed for the fitting of children under the age of 3 years or persons with a developmental age of under 3 years:

### **WARNING**

(for hearing aids NOT intended for small children and if conversion kit for child-safe battery door is available)

Choking hazard posed by small parts.

This hearing aid is not intended for the fitting of infants, children under 3 years and persons of mental incapacity.

### **WARNING**

for hearing aids for small children:

Instructions for the fitter of hearing aids for small children (e.g. in data sheet or fitting guide):

Choking hazard posed by small parts!

- ▶ Infants, children under 3 years and persons of mental incapacity must not wear the hearing aid without appropriate supervision.
- ▶ Acoustic coupling tube and earpiece must not be detachable by the child.
- ▶ The diameter of the earpiece is recommended to be in the range of 32 mm and may enclose the outer ear.
- ▶ Demonstrate and explain the battery door locking mechanism to the child's parents or caretaker.

# Important information

## Operating, transport and storage conditions

### For hearing aids with built-in power cells (lithium-ion rechargeable batteries)

Read the operating, transport and storage conditions in the hearing aid's user guide.

### For hearing aids with replaceable batteries

The hearing aids operate in the following environmental conditions (also valid in between uses):

|                      |                           |
|----------------------|---------------------------|
| Temperature          | 0 to 50 °C (32 to 122 °F) |
| Relative humidity    | 5 to 93 %                 |
| Atmospheric pressure | 700 to 1060 hPa           |

During extended periods of transport and storage, please observe the following conditions:

|                      | Storage                       | Transport                      |
|----------------------|-------------------------------|--------------------------------|
| Temperature          | 10 to 40 °C<br>(50 to 104 °F) | -20 to 60 °C<br>(-4 to 140 °F) |
| Relative humidity    | 10 to 80 %                    | 5 to 90 %                      |
| Atmospheric pressure | 700 to 1060 hPa               | 700 to 1060 hPa                |

For other parts, such as batteries, other conditions may apply.

## **Expected service lifetime**

The hearing aids are designed for an expected service life of 5 years. They are expected to remain safe for use for this duration.

## **Disposal information**

- ▶ Recycle hearing aids, accessories and packaging according to local regulations.

### **For hearing aids with replaceable batteries**

- ▶ To avoid environmental pollution, do not throw batteries into household trash.
- ▶ Recycle or dispose of batteries according to local regulations or return them to your Hearing Care Professional.

## **Conformance information**

The CE mark indicates conformity with the following European directives:

- Regulation (EU) 2017/745 for products manufactured by WSAUD A/S and 93/42/EEC for products manufactured by Signia GmbH and Sivantos GmbH\*
- 2011/65/EU RoHS concerning the restriction of hazardous substances
- Only for products with wireless functionality: 2014/53/EU RED concerning radio equipment

\* see declaration of conformity

The full text of the declaration of conformity can be obtained from [www.wsaud.com](http://www.wsaud.com).

The ACMA compliance mark  indicates conformity with the electromagnetic interference standards set by the Australian Communications and Media Authority (ACMA).

Devices with the FCC marking comply with the standards of the FCC regarding electromagnetic interference (only for products with wireless functionality).

### **Technical details**

If the certification information for USA is not listed in the hearing aids' user guide, the following certification numbers apply:

- Wireless platform, model: e2e 3.0
- **USA** FCC ID: SGI-WL003BTE

### **Notices**

This Class B digital apparatus complies with Canadian ICES-003.

Changes or modifications made to this equipment not expressly approved by the legal manufacturer may void the FCC authorization to operate this equipment.

This device complies with Part 15 of the FCC Rules and with ISED's licence-exempt RSSs.

Operation is subject to the following conditions:

- this device may not cause harmful interference, and
- this device must accept any interference received, including interference that may cause undesired operation.

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment to an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help.

For body worn operation, this device has been tested and meets the FCC RF exposure guidelines when used with the legal manufacturer's accessories supplied or designated for this product. Use of other accessories may not ensure compliance with FCC RF exposure guidelines.

Data is stored inside the hearing aid. This contains:

- the audiogram,
- the average daily wearing,
- the frequency distribution of acoustic environments (i.e., what percentage of the wearing time you were in a quiet environment, in a noisy environment, in an environment with speech, in a car and how often you listened to music),
- the different noise levels (including wind noise) within the last 24 hours,
- the average frequency of changes to user controls,
- the average frequency that algorithms in hearing aids are actively in operation,
- an identification number.

Part of the information inside your hearing aids is additionally stored in separate databases in anonymized form (meaning, it cannot be connected to you) in order to improve the quality of future hearing devices. This information may help your acoustician provide the optimal fitting for you. However, he/she will not read out this information without your consent. Also, when you send in your device for service, this information will not be read out without your consent.

## Nearfield magnetic induction functionality

The following table shows the main information on NFMi connectivity, some parameters are different from product to product, so the worst case is indicated:

|                         |                                     |
|-------------------------|-------------------------------------|
| Wireless technology     | Nearfield magnetic induction (NFMi) |
| Magnetic field strength | 0.14 A/m                            |
| Output power (EIRP)     | 53 $\mu$ W                          |
| Range                   | < 20 cm between hearing aids        |
| Center frequency        | 3.27 MHz                            |
| Channel                 | Single Channel Radio                |
| Bandwidth               | 330 kHz (-15 dB)                    |
| Modulation type         | PSK (Phase Shift Keying)            |

## Bluetooth® low energy\*

The following table shows the main information on BLE, some parameters are different from product to product, so the worst case is indicated:

|                     |                                              |
|---------------------|----------------------------------------------|
| Wireless technology | Radio frequency - Bluetooth low energy (BLE) |
| Output power (EIRP) | 3.9 mW                                       |

\* The Bluetooth word mark and logos are owned by the Bluetooth SIG, Inc., and any use of such marks by the legal manufacturer of this product is under licenses. Other trademarks and trade names are those of their respective owners.

|                 |                                                          |
|-----------------|----------------------------------------------------------|
| Range           | < 10 m between smartphone/<br>accessory and hearing aids |
| Frequency       | 2,400 – 2,480 MHz                                        |
| Channel         | 40 channels                                              |
| Bandwidth       | 2 MHz per channel ,<br>2 MHz spacing                     |
| Modulation type | Gaussian Frequency Shift Keying                          |

## EMI/EMC compliance

Wireless hearing aids comply with the following EMC/EMI standards:

| Standard                        | Test Type                                         | Note                                                                                                   |
|---------------------------------|---------------------------------------------------|--------------------------------------------------------------------------------------------------------|
| 47 CFR<br>Part 15,<br>Subpart C | RF emissions                                      | U.S. FCC requirements<br>for intentional radiators.                                                    |
| EN 300 330                      | RF emissions<br>including<br>spurious<br>emission | EMC and radio spectrum<br>matters for short range<br>devices in the frequency<br>range 9 kHz - 25 MHz. |
| EN 301 489-<br>1/3/17           | Immunity,<br>RF and ESD                           | Standard for low<br>power transmitters in<br>the frequency range<br>9 kHz - 40 GHz.                    |

| Standard       | Test Type                      | Note                                                                                                                                   |
|----------------|--------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
| EN 300 328     | Signal integrity and emissions | Wideband transmission systems; data transmission equipment operating in the 2.4 GHz ISM band and using wide band modulation techniques |
| IEC 60118-13   | RF immunity                    | International product standard for hearing aids to ensure adequate immunity to radio interference from mobile telephones.              |
| ANSI C63.19    | RF immunity                    | American National Standard method of measurement of compatibility between wireless communication devices and hearing aids.             |
| ANSI/AAMI PC69 | RF emissions                   | Implantable medical device EMC immunity, American National Standard.                                                                   |
| ISO 14117      | RF emissions                   | Implantable medical device EMC immunity, International Standard.                                                                       |
| EN 45502-2-1   | RF emissions                   | Particular requirements for pacemakers.                                                                                                |

|                                             |                                                                                                                           |
|---------------------------------------------|---------------------------------------------------------------------------------------------------------------------------|
| Emissions                                   | CISPR 11, Group 1, Class B                                                                                                |
| Immunity                                    | ESD, Conducted up to +/-8 kV,<br>Air discharge up to +/-15 kV<br><br>Power frequency magnetic fields:<br>30 A/m, 50/60 Hz |
| Immunity to RF fields                       | 80 MHz-2.7 GHz 10 V/m,<br>2.7 GHz-6 GHz 3 V/m                                                                             |
| Proximity to wireless communication devices | Up to 28 V/m –<br>See IEC 60601-1-2:2014/AMD1:2020,<br>Table 9                                                            |
| Proximity to magnetic fields                | 8 A/m@30 kHz, 65 A/m@134.2 kHz,<br>7.5 A/m@50 kHz                                                                         |

## **For US and Canadian markets only**

### **For prospective hearing aid users**

Good health practice requires that a person with a hearing loss have a medical evaluation by a licensed physician (preferably a physician who specializes in diseases of the ear) before purchasing a hearing aid. Licensed physicians who specialize in diseases of the ear are often referred to as otolaryngologists, otologists or otorhinolaryngologists. The purpose of medical evaluation is to assure that all medically treatable conditions that may affect hearing are identified and treated before the hearing aid is purchased.

Following the medical evaluation, the physician will give you a written statement that states that your hearing loss has been medically evaluated and that you may be considered a candidate for hearing aids. The physician will refer you to a Hearing Care Professional for a hearing aid evaluation.

The Hearing Care Professional will conduct a hearing aid evaluation to assess your ability to hear with and without hearing aids. The hearing aid evaluation will enable the Hearing Care Professional to select and fit hearing aids to your individual needs.

If you have reservations about your ability to adapt to amplification, you should inquire about the availability of a trial-rental or purchase-option program. Many Hearing

Care Professionals now offer programs that permit you to wear hearing aids for a period of time for a nominal fee after which you may decide if you want to purchase them.

U.S. federal law restricts the sale of hearing aids to those individuals who have obtained a medical evaluation from a licensed physician. U.S. federal law permits a fully informed adult to sign a waiver statement declining the medical evaluation for religious or personal beliefs that preclude consultation with a physician. The exercise of such a waiver is not in your best health interest and its use is strongly discouraged.

## **For hearing aid users**

Hearing aids will not restore normal hearing and will not prevent or improve a hearing impairment resulting from organic conditions. In most cases infrequent use of hearing aids prohibits the wearer from attaining the full benefit from it. The use of hearing aids is only part of hearing rehabilitation and may need to be supplemented by auditory training and instruction in lip reading.

### **Health considerations**

If soreness or skin irritation develops, discontinue wearing your hearing aids, and bring the hearing aids and earmolds to your Hearing Care Professional. Minor fit adjustments or earmold modification can often correct this condition. If soreness persists, discontinue wearing the hearing aid and see your physician. If excessive earwax accumulates when wearing your hearing aids, consult your Hearing Care Professional.

### **Battery tips**




If a battery is accidentally swallowed, seek medical attention immediately, or call the National Battery Ingestion Hotline collect at **1-800-498-8666**.

## **Identification information**

Your hearing aids have a serial number imprinted on them. The location of the serial number will vary according to the style of hearing aid you have chosen. Record the serial number in your user guide for future reference.


The year of manufacture is incorporated into the serial number.

For custom hearing aids, the year of manufacture is derived from the first two digits. In the following examples, the year of manufacture is 2022.

- ① 22: year manufactured
- ② X: facility code where hearing aid was manufactured
- ③ Serial number
- ④ L or R: left or right ear hearing aid
- ⑤ Manufacturer/Distributor
- ⑥ Model

For BTE and RIC hearing aids, the year of manufacture is derived from the second digit:

| Code (first digit) | Month | Code (first digit) | Month |
|--------------------|-------|--------------------|-------|
| A, N               | Jan   | G, U               | Jul   |
| B, P               | Feb   | H, V               | Aug   |
| C, Q               | Mar   | J, W               | Sep   |
| D, R               | Apr   | K, X               | Oct   |
| E, S               | May   | L, Y               | Nov   |
| F, T               | Jun   | M, Z               | Dec   |

| Code (second digit) | Year | Code (second digit) | Year |
|---------------------|------|---------------------|------|
| H, J                | 2025 | W, X                | 2020 |
| F, G                | 2024 | U, V                | 2019 |
| D, E                | 2023 | S, T                | 2018 |
| B, C                | 2022 | Q, R                | 2017 |
| Y, A                | 2021 | N, P                | 2016 |

Please note the warranty is based upon the date of purchase, not the date of manufacture.

Your Audiologist or Hearing Care Professional can answer any questions you may have about the identifying code on your hearing aid.

## **Tinnitus feature**

Your hearing aid may be equipped with a tinnitus feature. Many tinnitus patients also suffer from some degree of hearing loss. Many hearing aids can be used alone or in combination with the tinnitus feature. The tinnitus feature is fixed to a broadband noise that can be adjusted by your Hearing Care Professional for your tinnitus therapy.

The feature may provide temporary relief of your tinnitus.

## **Prescription Use Only**

U.S. federal law restricts this device to sale by or on the order of a doctor or Hearing Care Professional licensed to dispense hearing aids in your state. The use of any sound generating tinnitus therapy device should be only on the advice and in consultation with your Hearing Care Professional. Your Hearing Care Professional will properly diagnose and fit the device to your personal needs and requirements. This should include its use in a prescribed tinnitus treatment program. Your Hearing Care Professional will also be able to offer the appropriate follow-up care. It is important that you follow your Hearing Care Professional's advice and direction regarding such care.

There are some potential concerns associated with the use of any sound generating tinnitus therapy device. Discontinue use and seek medical evaluation if any of the following conditions occur:

- Chronic skin irritation on, near, or around the site of device placement.

- Unusual side effects (e.g., dizziness, nausea, headaches, heart palpitations).
- Perceived decrease in auditory function (e.g., decreased loudness, speech not as clear).

## **For Hearing Care Professionals**

### **Indications for use**

The tinnitus feature is a tool to generate sounds to be used in a tinnitus management program to relieve patients suffering from tinnitus. The target population is primarily the adult population over 21 years of age. The tinnitus feature is targeted for healthcare professionals treating patients suffering from tinnitus, as well as conventional hearing disorders. The fitting of the tinnitus feature should be done by a Hearing Care Professional participating in a tinnitus management program.

### **Device description**

The tinnitus feature is a software function that generates sound which is programmed into a hearing aid. Depending on the type of hearing aid, the hearing aid may be used in up to three modes of operation: as a hearing aid, as a tinnitus treatment device, or as a hearing aid and tinnitus treatment device.

When enabled, the tinnitus feature generates the sound and allows a patient's Hearing Care Professional to design and program appropriate settings for an individually prescribed sound treatment plan. The treatment plan should be used in a tinnitus management program for relief of tinnitus.

The tinnitus feature generates a broadband noise signal that varies in frequency and amplitude. These characteristics are adjustable by the Hearing Care Professional and are specific to the prescribed therapy designed by the professional for the patient's needs and comfort. The patient may have some control of the level or volume of the signal and the patient should discuss this adjustment as well as his or her comfort level and sound of the signal with their Hearing Care Professional.

A Hearing Care Professional should advise a prospective tinnitus feature user to consult promptly with a licensed physician (preferably an ear specialist) before using the tinnitus feature if the Hearing Care Professional determines through inquiry, actual observation, or review or any other available information concerning the prospective user that the prospective user has any of the following conditions:

- Visible congenital or traumatic deformity of the ear.
- History of active drainage from the ear within the previous 90 days.
- History of sudden or rapidly progressive hearing loss within the previous 90 days.
- Acute or chronic dizziness.
- Unilateral hearing loss of sudden or recent onset within the previous 90 days.

## **Cybersecurity**

Proper set-up and operation of your hearing aids requires a safe, private and secure operating environment.

The U.S. Federal Trade Commission's guidance for online security

<https://consumer.ftc.gov/online-security> provides examples of best practices for securing your home network and mobile device to ensure your security.

## **For children with hearing loss**

In addition to seeing a physician for a medical evaluation, a child with a hearing loss should be directed to an audiologist for evaluation and rehabilitation since hearing loss may cause problems in language development and the educational and social growth of a child. An audiologist is qualified by training and experience to assist in the evaluation and rehabilitation of a child with a hearing loss.

## **Warranty and service**

Your hearing aid, with the exception of the battery, is covered by a comprehensive warranty. All covered hearing aid parts received for warranty service at an authorized service center will be repaired or replaced with new or reconditioned components, without charge, to meet the performance specifications for that model.

This warranty does not cover malfunctions due to unusual wear and tear or mistreatment of the hearing aid such as physical shock, excessive wax build-up, or tampering with the hearing aid, any of which voids all warranties. Your Hearing Care Professional may charge a service fee for processing warranty service.

Warranty service must only be performed by an authorized service center. Service performed by unauthorized service depots voids this warranty, and repairs so necessitated will be done on a parts and labor cost basis.

Please refer to the warranty card included with your hearing aids for warranty period effective dates.

Your hearing aids may have additional loss and damage coverage. Please consult your Hearing Care Professional to determine if this is applicable to your hearing aids.

## **Support**

For more detailed information or in case of unexpected operation or events please contact

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© 09.2023, WSAUD A/S. All rights reserved
Sivantos GmbH and by Signia GmbH.
Sivantos GmbH y por Signia GmbH.