


Mojo MINI RIC RECHARGEABLE (V.2 RIC P-RECEIVER)

Moderate to severe-to-profound hearing losses
M7/M6/M5/M4/M3



Fitting range



Key Features	M7	M6	M5	M4	M3
Processing and fine-tuning channels	15	12	10	6	6
In-Situ Sound Optimisation	•	•	•	•	•
Experience Level Manager	•	•	•	•	•
Programs	5	4	3	3	3
Automatic Adaptive Directionality channels	15	12	10	6	6
Volume Control II*	•	•	•	•	•
Programmable Push Button**	•	•	•	•	•
Ambient Noise Rejection II	•	•	•	•	•
Noise reduction	•	•	•	•	•
Preference Log	•	•	•	•	•
Linked ZEN	•	•	•	•	•
Linear Frequency Transposition	•	•	•	•	•
Acoustic Situation Technology	11 (Linked)	7 (Linked)	4	3	
Impulse Noise Rejection II	•	•	•		
Automatic Adaptive Directionality	•	•	•		
Speech Enhancer II	II/Linked	Linked			
Real Ear Effect	•	•			
Wind Noise Control II	•				
High-frequency boost	•				
Connectivity					
2.4 GHz control (Android and iOS)	•	•	•	•	•
2.4 GHz audio streaming (iOS)	•	•	•	•	•
DEX assistive listening device compatibility	•	•	•	•	•
Accessories					
Mojo App	•	•	•	•	•
Acoustic remote control App	•	•	•	•	•
COM-DEX App	•	•	•	•	•

*Via DEX and app **Programmable: Volume control II, program shift or a combination of the two


Colours

114/Honey Blonde, 115/Chestnut Brown, 118/Tech Black, 119/Silver Grey, 120/Golden Brown, 121/Titanium Grey, 123/Autumn Beige


TECHNICAL DATA		EAR SIMULATOR IEC 60118-0:1983, + A1:1994	2CC COUPLER IEC 60118-0:2015, ANSI S3.22-2014
OSPL90	1600 Hz Peak Average	131 dB SPL 133 dB SPL 127 dB SPL	122 dB SPL 123 dB SPL 121 dB SPL
Acoustic output (Input 60 dB SPL)	1600 Hz Peak Average	116 dB SPL 119 dB SPL 108 dB SPL	107 dB SPL 109 dB SPL 104 dB SPL
Full-on gain (Input 50 dB SPL)	1600 Hz Peak Average	68 dB SPL 71 dB SPL 69 dB SPL	59 dB SPL 65 dB SPL 60 dB SPL
Acoustic Frequency range		100 Hz - 8050 Hz	100 Hz - 7550 Hz
Harmonic distortion (Input 70 dB SPL)	500 Hz 800 Hz 1600 Hz	<2% <2% <2%	<2% <2% <2%
Equivalent input noise		17 dB SPL	20 dB SPL
Battery current (stand by)		1.27 mA	1.27 mA
Battery current		1.43 mA	1.45 mA
Battery life / hours (Type 312 Zn-Air, 170 mAh)*		35 h	35 h
Mobile phone immunity (IEC 60118-13:2016, ANSI C63.19:2011)		IRIL: -25/-28/-10 dB SPL	U-rating: M4

*Battery life in real-life situations depends among other things on the hearing aid features used, streaming time, and the quality of the battery used.


MAXIMUM OUTPUT - EAR SIMULATOR



MAXIMUM OUTPUT - 2CC COUPLER



OUTPUT - EAR SIMULATOR



OUTPUT - 2CC COUPLER

