


Mojo CIC-M TR

Minimal to moderately hearing losses
M7/M6/M5/M4/M3/M2/M1



Fitting range



| Key Features | M7 | M6 | M5 | M4 | M3 | M2 | M1 |
|-------------------------------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
| Processing and fine-tuning channels | 15 | 12 | 10 | 6 | 6 | 4 | 4 |
| In-situ Sound Optimisation | • | • | • | • | • | • | • |
| Experience Level Manager | • | • | • | • | • | • | Manual |
| Programs | 5 | 4 | 3 | 3 | 3 | 3 | 1 |
| Ambient Noise Rejection II | • | • | • | • | • | • | • |
| Noise reduction | • | • | • | • | • | Minimal | Minimal |
| Preference Log | • | • | • | • | • | • | • |
| Volume Control II* | • | • | • | • | • | • | • |
| Zen | • | • | • | • | • | • | • |
| Linear Frequency Transposition | • | • | • | • | • | • | • |
| Acoustic Situation Technology | 11 | 7 | 4 | 3 | • | • | • |
| Impulse Noise Rejection II | • | • | • | • | • | • | • |
| Speech Enhancer II | II | • | • | • | • | • | • |
| High-frequency boost | • | • | • | • | • | • | • |
| Accessories | | | | | | | |
| Acoustic remote control App | • | • | • | • | • | • | • |

*Via app

Colours


025/Light Beige, 028/Clay Brown

TECHNICAL DATA


| | | EAR SIMULATOR IEC 60118-0:1983, + A1:1994 | 2CC COUPLER IEC 60118-0:2015, ANSI S3.22-2014 |
|---|-----------------------------|--|--|
| OSPL90 | 1600 Hz Peak Average | 109 dB SPL 118 dB SPL 109 dB SPL | 100 dB SPL 107 dB SPL 101 dB SPL |
| Acoustic output (Input 60 dB SPL) | 1600 Hz Peak Average | 94 dB SPL 103 dB SPL 90 dB SPL | 84 dB SPL 91 dB SPL 84 dB SPL |
| Full-on gain (Input 50 dB SPL) | 1600 Hz Peak Average | 58 dB SPL 60 dB SPL 57 dB SPL | 49 dB SPL 53 dB SPL 49 dB SPL |
| Acoustic Frequency range | | 100 Hz - 10000 Hz | 100 Hz - 8600 Hz |
| Harmonic distortion (Input 70 dB SPL) | 500 Hz 800 Hz 1600 Hz | <2% <2% <2% | <2% <2% <2% |
| Equivalent input noise | | 26 dB SPL | 26 dB SPL |
| Battery current (stand by) | | 0.91 mA | 0.91 mA |
| Battery current | | 0.93 mA | 0.95 mA |
| Battery life / hours (Type 10 Zn-Air, 100 mAh)* | | 110 h | 105 h |
| Mobile phone immunity (IEC 60118-13:2016, ANSI C63.19:2011) | | IRIL: -38/-32/-27 dB SPL | U-rating: M4 |

*Battery life in real-life situations depends among other things on the hearing aid features used, streaming time, and the quality of the battery used.


MAXIMUM OUTPUT - EAR SIMULATOR



MAXIMUM OUTPUT - 2CC COUPLER



OUTPUT - EAR SIMULATOR



OUTPUT - 2CC COUPLER

